EN S-CHECK
In vitro test for differentiation of vaccine strains from wild-type strains of *Salmonella* spp.

GENERAL INFORMATION

S-Check consists of two special media (medium A and B) for differentiation of *Salmonella* wild-type strains from the following vaccine strains of the same serovar:

- *Salmonella* Typhimurium mutant, strain 421/125, genetically stable, double-attenuated (histidine-adenine auxotrophic)
- *Salmonella* Enteritidis mutant, strain 441/014, genetically stable, double-attenuated (histidine-adenine auxotrophic)

BACTERIOLOGY

Generally, wild-type strains of the serovars Typhimurium and Enteritidis from *Salmonella enterica* subsp. *enterica* show a clear growth (turbidity) in both media A and B of S-Check. The vaccine strains only grow in medium B but not in medium A because of their special requirements for nutrients.

In very few cases the appearance of auxotrophic wild-type strains is observed. These wild-type strains do not grow either in medium A or in medium B.

INDICATIONS

The in-vitro test is indicated for use where it is necessary for epidemiological reasons to distinguish wild-type strains of *Salmonella* Typhimurium and *Salmonella* Enteritidis from the corresponding vaccine strain.

COMPOSITION
Medium A

Various salts, Glucose and Nicotinamide

Medium B

Various salts, Glucose, Nicotinamide and nutrients depending on the vaccine strain.

USE

After cultivation of the *Salmonella* spp. according to specified investigation method, selected colonies of typical morphology are transferred from the selective agar to Standard I Nutrient Agar or Columbia Blood Agar. These pure cultures are the basis for the differentiation with S-Check.

A small amount of colony material (1/4 or 1/2 colony, depending on size) is transferred via an inoculating loop (needle) to medium A and medium B. After transfer of the colony material, neither medium show any visible turbidity (see Figure 1).

Both media are incubated for 48 hours at 37°C.

To avoid contamination of the medium, please note the following:

- Do not use colonies taken directly from the selective growth medium. Always use pure cultures cultivated on Standard I Nutrient Agar or Columbia Blood Agar.
- Always use a laminar flow cabinet when transferring culture material.

After opening the bottles with medium, place the stoppers on its top. Do not touch the bottom part of the stoppers with your hands.

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COMPOSICIÓN
Medio A

Sales varias, Glucosa, and Nicotinamida

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
Evaluation (see also Figure 2 and 3):

Strain	Medium A	Medium B
Auxotrophic wild-type strain	No growth	No growth
Vaccine strain	No growth	Growth
Wild-type strain	Growth	Growth

EN : Figure 1: Medium A and B after inoculation (see arrow). You should be able to read a text when looking through.

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NOTE

1. S-Check may be left in their original bottles for test purposes. Following aseptic inoculation, the bottles may be sealed with their original stoppers.
2. After the test all the components of S-Check must be disposed according to professional standards.
3. Possible errors during use:
 - Failure to use a pure culture
 - Medium A and / or medium B show a distinct turbidity after inoculation (before incubation) – too high concentration of bacteria.
4. In the very few cases of finding auxotrophic wild-type strains the test should be repeated.
5. S-Check is not suitable for differentiation of mucoid *Salmonella* spp. colonies that may in very few cases be found on agar plates used for cultivation. Please contact the technical service or use other suitable methods (real time PCR, PFGE) for differentiation.

STORAGE

Store protected from light between + 2 °C and + 8 °C.

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
EN : Figure 2: Vaccine strain (see arrow). No growth in Medium A, growth in Medium B.

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
EN : Figure 3: Wild-type strain (see arrow). Growth in Medium A, growth in Medium B.

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En Colombia: Ceva Salud Animal S.A.S, Bogotá. Reg. ICA xxxx